

# The broadest collection of QC microorganism products for Growth Promotion Testing!



	<b>EZ-Accu Shot™ Select</b>	<b>EZ-Accu Shot™</b>	<b>EZ-CFU™ One Step</b>	<b>EZ-CFU™</b>
<b>AVAILABLE STRAINS</b>	5 compendial strains in 1 convenient kit	23 licensed, traceable strains	33 licensed, traceable strains	29 licensed, traceable strains
<b>NUMBER OF TESTS</b>	<ul style="list-style-type: none"> <li>• 10 tests per vial</li> <li>• 10 tests for each organism per kit</li> </ul>	<ul style="list-style-type: none"> <li>• 10 tests per vial</li> <li>• 50 tests per kit</li> </ul>	<ul style="list-style-type: none"> <li>• 19 tests per vial</li> <li>• 190 tests per kit</li> </ul>	<ul style="list-style-type: none"> <li>• 90+ tests per vial</li> <li>• 900+ tests per kit</li> </ul>
<b>HYDRATION FLUID (included in kit)</b>	Equilibrate to room temp	Equilibrate to room temp	Pre-warm hydrating fluid 30 minutes	Pre-warm hydrating fluid 30 minutes
<b>LYOPHILIZED MICROORGANISM PELLETS</b>	Use 1 pellet per suspension to obtain 10-100 CFU / 0.1 ml	Use 1 pellet per suspension to obtain 10-100 CFU / 0.1 ml	Use 2 pellets per suspension to obtain 10-100 CFU / 0.1 ml	Use 2 pellets per suspension to obtain 10-100 CFU / 0.1 ml
<b>RESUSCITATION OF MICROORGANISM</b>	Quick-dissolve pellet, no pre-incubation	Quick-dissolve pellet, no pre-incubation	Incubate suspension for 30 minutes before inoculation	Incubate suspension for 30 minutes before inoculation
<b>DILUTIONS</b>	No dilution step	No dilution step	No dilution step	One log dilution required
<b>STABILITY OF HYDRATED SUSPENSION</b>	8 hour stability*	8 hour stability*	8 hour stability*	30 minute stability

Choose from 4 unique products that meet your specific needs – all designed to deliver 10-100 CFU as required by the USP/Ph. Eur./JP.

[www.microbiologics.com](http://www.microbiologics.com)

Contact us to learn more  
320.253.1640  
800.599.BUGS(2847)

 **Microbiologics®**  
A safer, healthier world.