

# SterilAmp® II 5230

## Mesa Biological Indicators

Liquid Submersible SCBI for Low Temperature Steam

SterilAmp® II "5230" is a compact, self-contained ampoule biological indicator for monitoring industrial low temperature (110 - 118 °C) steam sterilization of liquids in small containers.

- **72-hour results**
- Small, conveniently-shaped ampoule
- No laboratory required
- Dramatic color change
- Packaged with 10 negative controls

26mm x 6.5mm



PASS

FAIL

SterilAmp® II "5230" contains a specially formulated culture medium that turns a dramatic yellow when spores grow, providing easy visual confirmation of sterilization within 72 hours, without the need for difficult and time-consuming laboratory work.

### Storage

Refrigerate at 2-8 °C.

### Certification

SterilAmp® II 5230 biological indicators are traceable to a recognized culture collection and certified for population, D-value, z-value, Fo value, and performance.

## Ordering Information

Species	Catalog #	Min. Population	Qty/Box
<i>B. subtilis</i> "5230"	SASU/6	1.0 x 10 <sup>6</sup>	100 Tests / 10 NC